

Tibial Tuberosity Advancement (TTA)



The **XORPTION** line of cages, used in **Tibial Tuberosity Advancements**, are manufactured using a specially selected compound of proven bioabsorbable materials. Material selection for each of our bioabsorbable implants is the most critical part of the product design. Each implant is evaluated independently with special consideration given to where it will be installed, the type and amount of loads placed on the implant, the healing profile for the fixation and the number of cycles placed on the implant. To ensure the proper material selection our engineers need to fully evaluate the indication and all the above parameters. The **XORPTION** Cage shares many of the same benefits of the **AdapTTA** cage.

XORPTION

Advantages of the XORPTION Cage and Quick Clip:

- Made from proven high-quality bioabsorbable polymers
- Engineered to provide adequate mechanical properties for safe bone healing
- Bioabsorbable profile mimics that of natural bone healing
- No tissue reaction & easy to remove if necessary
- Proprietary compounding of bioabsorbable materials with osteo-conductive materials provides faster bone growth
- Great scaffold for tissue growth
- Radiolucent, better postoperative imaging
- Available in 6mm, 9mm, 10.5mm, and 12mm sizes with Interchangeable wings



For more information, please contact us:

Everost US 855.EVEROST(383 7678)
Everost UK +44 1323 416057
info@everost.com
www.everost.co.uk

National Veterinary Services
T +44 1782 771100
E salesandmarketing@nvs-ltd.co.uk
W www.vetwholesaler.co.uk



EVEROST
advancing science, safety and innovation



National Veterinary Services

Tibial Tuberosity Advancement (TTA)



The **AdapTTA** line of cages, used in **Tibial Tuberosity Advancements**, are manufactured using the highest grade of proven biocompatible metals.

The innovative design provides adequate mechanical properties for safe and effective bone healing. With the novel "quick clip" the surgeon is no longer constrained to a specific wing size. The **AdapTTA** Cages allows the surgeon to pick the correct wing size for the patient and attach the wing onto the desired cage size.

Advantages of the AdapTTA Wing

- Interchangeable wings; choice made in surgery!
- 360 degree of wing freedom
- Available in biocomposites, titanium and stainless steel
- Limited contact / less implant material
- Available in 8 sizes
- Proprietary surface finish



CAGE & WING CHART

| CAGE & WING CHART | | | | | | |
|-------------------|----------|-----------|-----------|-----------|-----------|------------------|
| CAGE WIDTH | | | | | | |
| 16 | | | | | | X |
| 12 | | | | X | | X |
| 10.5 | | | X | X | | X |
| 9 | | X | X | X | | X |
| 7.5 | X | X | X | X | | X |
| 6 | X | X | X | X | X | X |
| 4.5 | X | X | X | X | X | X |
| 3 | X | X | X | X | X | X |
| | 9 | 11 | 13 | 15 | 17 | 19 |
| | | | | | | WING SPAN |

For more information, please contact us:

Everost US 855.EVEROST(383 7678)
 Everost UK +44 1323 416057
 info@everost.com
 www.everost.co.uk

National Veterinary Services
 T +44 1782 771100
 E salesandmarketing@nvs-ltd.co.uk
 W www.vetwholesaler.co.uk



EVEROST
 advancing science, safety and innovation

